

Number order confirmation:

# What to consider before setup?

In order to set up the device correctly and in time, it is important to inspect the on-site installation requirements in advance using this check list. If necessary consult an expert.

If the installation is done by the MELAG Customer Service, send this list to MELAG (email: info@melag.de or fax +49 (0)30 757 911 1083).

Practice / clinic (name, address / stamp)	Stockist (name, address / stamp)

X Tick the appropriate items on selection or select the applicable one:

 $\checkmark$ 

present/fulfilled



not present/not fulfilled

### Space requirements / premises

Note the following:

- ▶ If the device is installed directly in the corner of an L-shaped work surface, existing drawers to the left or right of the device **cannot** be opened.
- If the device is installed with a floor unit, it must be additionally secured to the wall using the wall fixing (included in the scope of delivery of the floor unit).

	✓ ×
The width of the door etc. from the practice entrance up to the setup site is at least 60 cm.	$\circ$
The space required for MELAtherm is available.	$\circ$
The additional space required for a mixed-bed resin cartridge is available.	$\circ$
The installation surface can withstand a floor load of 113 kg for an under-desk unit, 119 kg for a free-standing device and 182 kg for a top-frame device.	$\circ$

Space requirement MELAtherm 10/MELAtherm 10 Evolution			
Width x Depth	60 x 60 cm		
Height	min. 82 cm		
Height with stainless steel cover plate	84 cm		
Height with floor unit	124 cm		
Distance to opposite work surface	min. 70 cm		

Additional space required for a mixed-bed resin cartridge	MELAdem 53	MELAdem 53 C
Diameter x Height	`	24 x 45 cm (approx. 49 cm incl. pipe elbow)

		lis	

MELAtherm 10/MELAtherm 10 Evolution

Number order confirmation:

# **MELAG**

#### **Electrical connection**

4	Δ	
	Ĭ	/

## **WARNING**

Danger of electric shock

Work on the electrical system should only be performed from an authorised and qualified electrician.

			×			
Th	The practice operator provides a correctly installed electrical installation.  The connection should be made with cable cross-sections of 2.5 mm².					
	MELAtherm 10 DTA/MELAtherm 10 Evolution DTA:  CEE socket with 380-415 V 50 Hz or 60 Hz*) 16 A fuse protection type B					
	E socket with 380-415 V, 50 Hz or 60 Hz*), 16 A fuse protection type B	$\bigcirc$	$\bigcirc$			
All	-pole disconnection is guaranteed	$\bigcirc$	$\bigcirc$			
Ad	Additional RCD switch with 30 mA tripping current is available					
ME	ELAtherm 10 DTB/MELAtherm 10 Evolution DTB:					
22	0-240 V, 50 Hz or 60 Hz*) socket, 16 A fuse protection type B	$\bigcirc$	$\bigcirc$			
Additional RCD switch with 30 mA tripping current is available						
*) ]	The 60 Hz device version only applies to MELAtherm 10, but not to MELAtherm 10 Evolution.					
	PLEASE NOTE  If the DTB device version is intended for permanent operation with an ambient temperature of 35 °C, then contact the technical service for further information.	over	,			
W	ater connection					
No	te the following:					
	The quality of the water used for reprocessing influences the value-retention of the load. MELAG recomme a pre-filter with a mesh size of ≤ 100 µm to prevent particles from entering from the water pipe (e.g. sand, I					
	If the water hardness is higher than the maximum permissible hardness, a water softening system must be upstream. When using a central demineralisation system, the maximum permissible conductivity is $5 \mu\text{S/cr}$		alled			
<b>•</b>	MELAG recommends a leakage water detector, e.g. the water stop from MELAG (art. no. ME01056).					
	Observe the country-specific regulations on drinking water protection, e.g. SVGW, KIWA, WRAS.					
		$\checkmark$	X			
	e water outlet was properly installed under observance of the technical manual. The water outlet is able to charge without any impediment (observe any pumping stations).	$\bigcirc$	0			
CA	UTION! The wastewater can reach up to 95 °C for short terms.					
Th	e pump-off rate amounts to 40 l/min, the amount of wastewater up to 30 l/h (in smaller intervals).	$\bigcirc$	$\bigcirc$			
Th	e maximum permissible water hardness amounts to 40 °dH.	$\tilde{\bigcirc}$	$\tilde{\bigcirc}$			
	nnections (3/4" external thread) to the drinking water are present and accessible for the connection of the ter treatment unit.	Ŏ	Ŏ			
	e recommended flow pressure for the cold water connection and the demineralised water connection is i-10 bar in the Netherlands).	$\bigcirc$	$\bigcirc$			

Pro	cess	age	ents
	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	June	,,,,

		$\checkmark$	×
	ice operator provides the process agents required for commissioning in accordance with the ecommendation (one 5 I container each of cleaner and neutraliser and at least one 1 I container of	0	0
Remark	rs and notes		
The follow	ing points are still open and must be completed by the installation date:		
	PLEASE NOTE All non present/not fulfilled points are to be fulfilled by the practice operator until installation points remain open when the device is installed and the installation cannot be carried out, the operator bears the costs for the service call.		lf
Other rema	arks:		
The inform	potion outlined above has been complied. The requirements for the installation and commissioning o	ftho	
	nation outlined above has been complied. The requirements for the installation and commissioning c guaranteed.	i the	

Signature / stamp practice / clinic

Signature / stamp customer service / stockist