

Please give the postal address of the specialist dealer/depot:

Specialist dealer/depot

Street/house number

Post code/town

Country

Dear doctor,

for validating your sealing processes MELAG provides a seal seam stability test for a cost price of 125.00 € plus V.A.T. (state 01/2015).

Please complete this form and send it including all film test strips in an A4 envelope with address window to your specialist dealer/depot.

After the film test strips have been tested, MELAG will issue verification of testing and if successful, a certificate. This certificate confirms conformity of the seal seam with the standard DIN EN 868-5, Appendix D.

Contact data

| Specialist dealer/depot | Practice/clinic |
|-------------------------------|-------------------------------|
| Specialist dealer/depot | Practice/clinic |
| Contact person | Contact person |
| Street/house number | Street/house number |
| Post code/town | Post code/town |
| Country | Country |
| E-Mail address | E-Mail address |
| Telephone for further queries | Telephone for further queries |

Technical Data


| Device description | Settings |
|--|--|
| <input type="checkbox"/> MELAseal [®] 100+ | |
| <input type="checkbox"/> MELAseal [®] 200 | Sealing temperature/thermostat setting |
| <input type="checkbox"/> MELAseal Pro [®] | |
| <input type="checkbox"/> other MELAG sealing device: | Sealing force in N ¹⁾ |
| MELAseal | Sterilized with (device description, manufacturer) |
| Serial number of sealing device | Serial number of the steam sterilizer |

¹⁾ not required if the sealing device is not validated and serial number and device description are filled in

Please turn over →

Please observe the following sequence:

1. Cut 3 x 10 cm length test strips from the smallest packaging format (e.g. MELAfo[®] 502) used in your practice.
2. Cut a 1 x 10 cm long test strip from every other larger format.
3. Please mark each test strip with a sequential number.
4. All test strips must have two seal seams separated by a gap of 5 cm.



5. Sterilize all test strips in your autoclave using the Universal-Program (134°C).
6. Please fill out the following table.

| Film no. | Manufacturer/designation | Width of test strip [in cm] | Use by/shelf-life date of the film ²⁾ |
|----------|--------------------------|-----------------------------|--|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

²⁾ You can find the date on the transport packaging.

Should you have any further questions regarding performing the seal seam stability test, please do not hesitate to contact us:

MELAG Medizintechnik oHG

Geneststraße 6-10
10829 Berlin
Germany

E-mail: info@melag.de
Web: www.melag.de

_____ Date

_____ Signature Client